Upoznajte METTLER TOLEDO

CROLAB sekcija za hranu

Koprivnica, 25.09. 2012

METTLER TOLEDO

Jeste li znali ...

... da vaša zubna pasta u prosjeku sadrži više od 50 sastojaka?

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Kalibracijski laboratorij / Calibration laboratory

Vage / Weighing instruments

Akreditiran od strane / Accredited by the

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This calibration certificate documents the traceability to national standards, which realize the units of measurement according to the International System of Units (SI).
The DKD is signatory to the multilateral agreements of the European co-operation for Accreditation (EA) and of the International Laboratory Accreditation Cooperation (ILAC) for the mutual recognition of calibration certificates.
The user is obliged to have the object recalibrated at appropriate intervals.

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This calibration certificate may not be reproduced other than in full except with the permission of both the Accreditation Body of the DKD and the issuing laboratory. Calibration certificates without signature and seal are not valid.

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E2

Kalibracija s normalnom rezolucijom
Calibration with normal resolution

Instrument je podešen prije kalibracije
The instrument was adjusted before calibration.

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
Kalibracija je provedena u skladu s OIML R111:2004 i OIML R111:2004. The calibration was performed according to OIML R111:2004 and OIML R111:2004.

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The net display values shown in the table are the minimum weights. At these values the measurement uncertainty of the display is equal to or less than 1/1 (no safety factor), 1/3, 1/5 or 1/1z of the required tolerance in the table.

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The measurement uncertainty values used are taken from the DKD Calibration Certificate and based on the expansion factor $k = 2$ (95% probability of occurrence).

Kalkulacije minimalnih neto prikaza (minimalnih odvaga) za različite tolerancije odvaga uz različite faktore sigurnosti.

Sample calculations of minimum net display values (minimum weights) for various weighing tolerances at different safety factors.


Odnos između prikazane mjerne nesigurnosti i zahtijevane mjerne nesigurnosti Relationship between the measurement uncertainty of the display and the required tolerance						
Mjerna opseg 1 / Range 1	Tražene tolerancije / Required tolerance	1/1 (no safety factor)	1/2 (safety factor of 2)	1/3 (safety factor of 3)	1/5 (safety factor of 5)	1/10 (safety factor of 10)
	0,1 %	0,00822 g	0,01654 g	0,02497 g	0,04216 g	0,08717 g
0,2 %	0,00410 g	0,00822 g	0,01236 g	0,02074 g	0,04216 g	0,04216 g
0,5 %	0,00164 g	0,00327 g	0,00492 g	0,00822 g	0,01654 g	0,01654 g
1 %	0,00082 g	0,00164 g	0,00245 g	0,00410 g	0,00822 g	0,00822 g
2 %	0,00041 g	0,00082 g	0,00123 g	0,00204 g	0,00410 g	0,00410 g
5 %	0,00016 g	0,00033 g	0,00048 g	0,00082 g	0,00164 g	0,00164 g

Korisnik je odgovoran za namještanje i održavanje postava (setting) vage koje su bile korištene za određivanje minimalne odvage vage.

The user is responsible for maintaining the configuration (setting) of the scale/balance which was used when the minimum weight was determined.

Datum kalibracije / Date of calibration: 2008-09-12

Osoba za osobu / Person in charge

Ivan Maras

Pet proizvodnih procesa

- Laboratorijska oprema
- Procesna oprema
- Industrijska oprema
- Kontrola proizvoda
- Trgovačke vage

Mettler Toledo vam nudi

Proizvodnja / higijenski dizajn

- Platforme, vage i terminali izrađeni prema standardima za prehrambenu industriju
- do IP69k zaštite za pranje pod tlakom za korištenje u higijenskoj zoni

Završna kontrola

- Integralna Statistička Kontrola Kvalitete (SQC) 
- Kontrolne vage, Metal detektore, uređaje za kontrolu sa X-zrakama

Slijedivost

- kompletna logistička rješenja (CargoScanneri)

Validacija / kvalifikacija / servis

- osigurava ispravan rad Vaših sistema
- nadzor nad uređajima prema ISO17025 standardu
- Zakonska ovjera mjerila

Vage za higijensku zonu

- Vage, platformske, mostne ili sa ćelijama za vaganje spremnika u proizvodnji
- Terminali za sve proizvodne aplikacije (doziranje, recepture, kontrole,)
- Povezivost u jedinstvene sustave sa nadzorom
- Rad prema standardima GMP, GAMP, IFS, BRC, HACCAP,....

Kontrola težine proizvoda

- Dinamička kontrolna vaga za kontrolu pakovina
- Do 400 odvaga u minuti
- Proizvodi od 5 g do 100 kg
- Odbacivanje proizvoda koji su van postavljenih tolerancija
- Povratna veza na punilicu
- Povezivanje u mrežu na FreeWeighNet program
- Do IP69K zaštite



Kontrola proizvoda, metal detektori

- Kontrola kontaminacije proizvoda metalom
- Kontrola pakovina, sirovina ili proizvoda prije pakiranja
- Sustav se projektira prema proizvodnoj liniji i aplikaciji
- Prepoznavanje Feritnih, neferitnih i inox metala

Kontrola proizvoda sa X-zrakama

- Kontrola X-zrakama, zasnovana na razlici gustoće proizvoda i kontaminata
- Kontrola kontaminacija metalom, staklom, porculanom, tvrdom plastikom....
- Kontrola pakovina, sirovina ili proizvoda prije pakiranja
- Kontrola oštećenosti i nedostatka proizvoda
- Pronalaženje stakla u staklu

Tipovi uređaja

Staklo u staklu

Kontrola kutija i vrećica

Kontrola kapljovina

Kontrola transportnog pakiranja



Higijenski dizajn

Kontrola sirovina


Povezivanje u mrežu

- Povezivanje u mrežu statičkih i dinamičkih vaga, kao i ostalih mjernih uređaja
- FreWeighNet program


Što je FreeWeigh.Net?

Sistem za kontrolu punjenja i kontrolu kvalitete omogućuje:

- **Fleksibilni dizajn na principu client-server mreže**
- **Kompatibilnost sa SQL bazom podataka**
- **Ethernet mreža za rad**
- **Podrška za više korisnika**
- **Fleksibilno i intuitivno sučelje**
- **Robusnost**
- **Kontrola težina, volumena, pH, dimenzija, kontaminacije, atributa, ...**
- **Kreiranje raporta prema vlastitim zahtjevima**

FreeWeigh.Net veliki sistem

Hvala na pažnji

METTLER **TOLEDO**